


A brief history of ICP & ICP/MS Nebulizers




Presented by John A. Burgener


**Pneumatic nebulizers
have always been,
and are still
the main
sample introduction method
for ICP and ICP/MS.**

John A. Burgener

 *Burgener Research Inc.*



In the beginning, there was the “Adjustable Cross Flow”



*They were difficult to adjust, easily broken
and easily plugged from salting or particles.*


John A. Burgener

 Burgener Research Inc.

This was soon followed by the Meinhard Concentric, 1974



Meinhard Type T nebulizer
Serial # 4250
Made approx. 1978



*Simple, elegant design.
Similar to modern Concentrics.
But salted and plugged easily.*

John A. Burgener

 Burgener Research Inc.


1947: My father started Technical Service Laboratories

**1974: TSL received the second ICP made by JA
an Atom comp 950 - Arc/Spark instrument
with an ICP source replacing the arc spark box.**

**1977: I was asked by my father to help get the ICP running.
The objective was to run Whole Rock analysis.
Whole Rock solutions were 3% LiBO4
ALL nebulizers salted in hours or minutes.
The Adjustable Cross Flow was the first one we used.
Meinhard's concentric salted too fast - minutes.
The cross flow salted in hours,
and plugged too easily.**

Induction pulls liquid into gas stream.

Low pressure beside gas stream causes liquids to vaporize leaving dissolved salts behind, causing salt deposits that plug gas orifices.



Concentrics and Cross flows plug easily. To minimize plugging, Babington designs were developed.


Jan. 14, 1969

R. S. BABINGTON ETAL

3,421,692

METHOD OF ATOMIZING LIQUIDS IN A MONO-DISPERSED SPRAY

Filed Dec. 29, 1966




John A. Burgener


 Burgener Research Inc.

**High Solids resistant, but very inefficient.
Only a little part of the liquid passed near the orifice.**

1978 Babington V Groove by Suddendorf and Boyer



1980 Babington Pillar and Post by Garbarino and Taylor




**High Solids resistant, but very inefficient.
Only a little part of the liquid passed near the orifice.**

John A. Burgener

 Burgener Research Inc.

~ 1982 GMK Nebulizer Glass Babington V Groove




**Handled particles well, but
gas orifice plugged if running high salts.**

John A. Burgener

 Burgener Research Inc.

1982 - Meddings' MAK Nebulizer A Fixed Cross Flow.

Significant improvement over previous ones.




Handled high salts better than most.
More stable, no adjustments.
But still easy to plug with particles.

John A. Burgener

 Burgener Research Inc.

1983 Meinhard introduces the C Type Nebulizer



**Lapped edges reduced salting,
but still salted on our solutions.**

John A. Burgener

 Burgener Research Inc.

1983 Precision Glassblowing began making glassware and nebulizers




**Similar to Meinhard's A Type
and still salted on our solutions.**

John A. Burgener

 *Burgener Research Inc.*

~ 1983 Jarrell Ash (Thermo) Sapphire V Groove Nebulizer




**Handled particles well, but
Sapphire gas orifice would still salt rapidly.
And Sapphire bead would often fall out.**

John A. Burgener

BR Burgener Research Inc.

1984 Meinhard introduces the K Type Nebulizer




Recessed inner capillary significantly reduced salting, but still salted on our solutions.

John A. Burgener

 Burgener Research Inc.

~1984 Glass Expansion began making ICP glassware



John A. Burgener

 Burgener Research Inc.

1984-85 The Burgener - Legere Nebulizer

First commercial Teflon Nebulizer.
V Groove with no adjustable parts.




Handled high salts,
More stable, no adjustments,
Not plugged with particles.
Very difficult to make - about a week each.
RAN FOR MONTHS WITHOUT SALTING!

John A. Burgener

BR Burgener Research Inc.

1986 Fassel, Rice, Lawrence: MicroNebulizer for Direct Injection (US Patent # 4,575,609)




One of the first micro-nebulizers
could not handle salts
was not useful for most samples.

John A. Burgener

 Burgener Research Inc.

1986 Hildebrand Grid Nebulizer



Excellent Sensitivity
Difficult to clean
Long washout times

John A. Burgener

 Burgener Research Inc.

Late 1980s or early 1990s Perkin Elmer Gem Tip Cross flow




Excellent general purpose nebulizer for the analysis of strong mineral acids (including HF) and samples with less than 5% dissolved solids.
Up to 20% NaCl can be aspirated by the nebulizer for 1 hour without clogging.
Uses chemically resistant GemTips made of sapphire and ruby in a Ryton end cap for maximum resistance.

1988 CETAC Ultrasonic Nebulizers

USN 5000, 6000

Efficiency approaches 30%, Sensitivity improves ~10x
Droplet size < 5 microns.
Potentially heavy solvent load so Desolvation essential.
Peltier cooling essential.
Membrane separator available and usually essential.
Desolvation interferences occur (eg., As III vs. As IV)

Does not handle high solids well.
Long washout times required due to desolvators.



Evaluation of an Ultrasonic Nebulizer Using Perkin-Elmer Sequential ICP Instrumentation - Robert Thomas, Cindy Anderau - Atomic Spectroscopy, 10 (2), 71 (1989). This paper describes the capabilities of the first commercially available ultrasonic nebulizer, the CETAC USN-5000 with two of Perkin Elmer's ICP optical emission spectrometers, the Plasma 40 and Plasma II. Various criteria were evaluated including detection limit, precision and memory test. The conclusion of the work was that the early design of this product showed great promise but required some improvements to be accepted and used as a truly routine analytical tool.

Late 1980s: Cyclonic Chambers Arrive

The early chambers were variations on AA chambers.
Most were large straight tubes with baffles.
The Scott Double Pass became the most successful of these.
Most nebulizers were designed to fit 35 mm OD front ends of chambers.


In the late 1980s cyclonic chambers became popular.
introduced by Glass Expansion

With 2X the sensitivity, and 1/2 the washout time,
most labs switched to cyclonic chambers.

At first, only concentrics fit.

So nebulizers were re-designed to fit the GE standard 6mm
cyclonic chamber nebulizer fittings.

1987 Glass Expansion's first Nebulizer: The VeeSpray a ceramic V groove Nebulizer



Material: alumina ceramic
High TDS tolerance: typically up to 30%
High tolerance to particulates: typically up to 300 μ m
Must be pumped - does not self aspirate
Operates best between 1.5 and 2.5mL/min

John A. Burgener

 Burgener Research Inc.

1989 Glass Expansion's first Concentric The Conikal




**Machined inner capillary
instead of glass blown.**

John A. Burgener

 Burgener Research Inc.

1989 Noordermeer Glass V Groove

Sample




**Good to Excellent Sensitivity
Orifice can salt.**

John A. Burgener

 Burgener Research Inc.

1992 Glass Expansion began making the non salting Sea Spray




**Tip of inner capillary is recessed more than with Meinhard K
Nose body extended to create a self cleaning area inside the tip
Can handle 20% salts without blocking
Had improved sample line attachment for easy fitting
and minimal dead volume.**

John A. Burgener

 Burgener Research Inc.

1993 Modified Lichte Glass V Groove




Good to Excellent Sensitivity
Orifice can salt.
Overloads many ICP torches.

John A. Burgener

 Burgener Research Inc.

1993 BTF - The first Burgener Parallel Path




Good to Excellent Sensitivity
Non salting design:
 can run 100% saturated solutions.
Lasts for months without washing.
Difficult to make - about 2 days each.
Not as sensitive as a concentric.

John A. Burgener

 Burgener Research Inc.


The Main Burgener Parallel Path Nebulizers 1994 - 1995



The BTS 50



The Thermo Trace



The T2002

**Less Difficult to make - about 1 day each.
Eventually we reduced it to about 4 hours each.**

Mid 1990s: Perkin Elmer Miniature V Groove The Gem Cone

Last V groove design I am aware of



John A. Burgener

 Burgener Research Inc.

V groove designs still commercially available:

As of 2013


Glass Expansion VeeSpray




Perkin Elmer Gemcone



**Precision Glassblowing
Modified Lichte**



**Precision Glassblowing
Noordermeer**



They all suffer from salting, and are no longer significant in the market. They have generally been replaced by newer, better designs.

Cross Flow nebulizers also declined in use.

Partly in response to the cyclonic chambers.

They are not as good as concentrics in sensitivity and stability.

They are less easy to salt and plug as a concentric, but it is still easy.

But they are inexpensive to make.

So most manufacturers have at least one model available.

And some new ones have been introduced in the last 10 years:

**SavilleX: X-Flow
PFA Cross Flow Nebulizer**



**SCP: Ultem Mini Cross Flow
designed for cyclonic chambers**



John A. Burgener

 Burgener Research Inc.

The arrival of the ICP / MS

In the late 1980s and early 1990s,
ICP / MS became useful lab instruments.

Previously they were mainly used by
R&D institutions and universities.


With Agilent and Perkin Elmer / Sciex
ICP /MS became common instruments by the mid 1990s.

Nebulizers running 2 ml/min sent too much sample to the Mass Spec.

Concentrics were able to run much lower flow rates,
and still provide excellent stability and sensitivity.

The development of micro nebulizers began.

MEINHARD HEN: High Efficiency Nebulizer ~ 1993



	TR-30	HEN-170
Sample Capillary:	220-320 μm	70-110 μm
Gas orifice ID:	350-450 μm	150-200 μm
Pressure at 1 L/min:	25-40 psi	150-180 psi
Normal operating pressure:	15-25 psi	25-60 psi

Very Low Flow capable
Very easily plugged and salted

Burgener Nebulizers Showed the advantages of Teflon



**The Non Salting nature of Teflon Gas Orifices
was desirable in all designs.**

1997 Cetac Microconcentric Nebulizer

The MCN 100 was probably the first Teflon Concentric.
It was also the first non-Burgener Nebulizer
commercially available that was 100% Teflon.




MCN 100 and Aspire
50, 100, 200, or 400 μ L/min

John A. Burgener

 Burgener Research Inc.

1997 Direct Injection High Efficiency Nebulizer Meinhard's DIHEN



DIHEN: inserted directly into a demountable torch

Normal carrier flow: 0.2 – 0.4 L/min

Normal Operations: 30 – 50 psi, 5 – 200 μ L/min

Natural aspiration rate: 10 – 300 μ L/min

Minimizes speciation interferences

Easily introduce highly volatile solvents

Essentially 100% transport

Design and test pressure: 90, 120, 150, or 170 psi @ 1 L/min

John A. Burgener

 Burgener Research Inc.

A. Montaser, J. A. McLean, and Jerold M. Kacsir,
A Novel Direct Injection High Efficiency Nebulizer for Atomic Spectroscopy, 1997
Patent Applications filed December 28, 1998
U.S. Patent No. 6,166,379, December 26, 2000.

1997 Direct Injection High Efficiency Nebulizer Meinhard's DIHEN



**Good Concept,
But droplets do not have time to vaporize.
droplets arrive in plasma as droplets.**

**Other designs will provide greater sensitivity and stability
by NOT being close to the plasma.**


John A. Burgener

 *Burgener Research Inc.*

A. Montaser, J. A. McLean, and Jerold M. Kacsir,
A Novel Direct Injection High Efficiency Nebulizer for Atomic Spectroscopy, 1997
Patent Applications filed December 28, 1998
U.S. Patent No. 6,166,379, December 26, 2000.

Note: 1997 - J E Meinhard Produces the 50,000th Nebulizer!

1999 Elemental Sciences introduced their PFA Concentric Nebulizers



20, 50, 100, or 400 $\mu\text{L}/\text{min}$

John A. Burgener

 Burgener Research Inc.

1999 Burgener Introduces the Micro 1



**Designed to run low flows, 50 to 1,500 μ l/min
New method of construction - New body design.
Still based on the original Parallel Path method.**


John A. Burgener

 Burgener Research Inc.



What's New in Nebulizers

PART 2: Since the year 2000



2000 Burgener Micro 3



**Improved Body design over Micro 1,
still original Parallel Path Method
Still aiming at low flow
50 - 1500 ul/min range
Looks like a Mira Mist but very different.**

John A. Burgener

 Burgener Research Inc.

More changes in PFA Concentrics then in Glass Concentrics

Elemental Scientific: Many models including micro PFA nebulizers
Glass Expansion: Opal Mist Nebulizer
Cetac: Aspire, MCN 100 Micro Concentrics
Savillex: C Flow

**Many micro designs, and standard flows
But nothing essentially unique or exciting since 2000**

Models range from a few microliters/minute to about 1 ml/min

ESI produces many variations on PFA Neb.

PFA MicroFlow Nebulizers

PFA-20	20-35 $\mu\text{L}/\text{min}$
PFA-50	45-65 $\mu\text{L}/\text{min}$
PFA-75	65-96 $\mu\text{L}/\text{min}$
PFA-100	90-150 $\mu\text{L}/\text{min}$
PFA-200	150-250 $\mu\text{L}/\text{min}$
PFA-400	300-400 $\mu\text{L}/\text{min}$

PFA-ST Nebulizers

0.15 mm i.d.	20 $\mu\text{L}/\text{min}$
0.20 mm i.d.	50 $\mu\text{L}/\text{min}$
0.25 mm i.d.	100 $\mu\text{L}/\text{min}$
- pumped	0.1 to 1.0mL/min
0.30 mm i.d.	200 $\mu\text{L}/\text{min}$
0.50 mm i.d.	400 $\mu\text{L}/\text{min}$
- pumped	0.5 to 3.0mL/min
0.80 mm i.d.	700 $\mu\text{L}/\text{min}$

PFA-HS: High dissolved solids, 300 μm ID


PFA-LS: laser ablation applications, 20 $\mu\text{L}/\text{min}$
produces an aerosol at low gas flow rates:
0.2—0.3 L/min

PFA-LC: HPLC Interface for micro and standard flow
Low internal volume, 150 μm ID



PFA-Q: internal quartz capillary,
for applications where the analyte
has memory issues with plastic
such as Iodine

PolyPro-ST Nebulizers

0.15 mm i.d.	20 $\mu\text{L}/\text{min}$
0.20 mm i.d.	50 $\mu\text{L}/\text{min}$
0.25 mm i.d.	100 $\mu\text{L}/\text{min}$
- pumped	0.1 to 1.0mL/min
0.30 mm i.d.	200 $\mu\text{L}/\text{min}$
0.50 mm i.d.	400 $\mu\text{L}/\text{min}$
- pumped	0.5 to 3.0mL/min
0.80 mm i.d.	700 $\mu\text{L}/\text{min}$
- pumped	up to 5.0mL/min



2012 Glass Expansion changes Opal Mist




Better sample line attachment
More secure gas line attachment
Less contaminating gas line attachment

John A. Burgener

 Burgener Research Inc.

2005 Savillex C-Flow and X-Flow



John A. Burgener

 Burgener Research Inc.

Refinements in Glass Concentrics

Most manufacturers have continuously improved their designs

John A. Burgener

 Burgener Research Inc.

1980 & 90s Meinhard Concentric

Back end was much larger than inside central capillary.
Led to plugging inside capillary,
And poor sample line connections.



John A. Burgener

 Burgener Research Inc.

1980 & 90s Meinhard Concentric

Sample line was pressed into place
Held partly by suction of nebulizer
and partly by the capillary sticking in place
and partly by luck.



Gas lines were Tygon tubing pushed over the barb.
they regularly popped off if the pressure was 40+ psi

Modern Meinhard Concentric

Back end is only a bit larger than inside central capillary.
Minimizes plugging inside capillary.
New Sample line fittings:
Secure, easy to attach to LC or pump.




New Gas lines are connected
with PFA fittings.
Tight fit - no lines popping off.
Non contaminating.

John A. Burgener

 Burgener Research Inc.

MEINHARD®

Comparison of Trace Metals Leached from Nebulizer Gas Fitting



Conventional gas line connector



PFA FG⁺ quick connect gas line

1990s Glass Expansion

Sample line was held with Teflon fitting.



Gas lines were Tygon tubing pushed over the barb.
they regularly popped off if the pressure was 40+ psi

John A. Burgener

 Burgener Research Inc.

2012 Glass Expansion

**Machined inner capillary “Vitricone design”
Constant ID for whole length:
plugging mainly occurs at connection.**




**Gas lines are connected
with secure fittings with O ring seals.
Tight fit - no lines popping off.**

John A. Burgener

 Burgener Research Inc.

2012 Glass Expansion



Machined inner capillary “Vitricone design”

Constant ID for whole length:


plugging mainly occurs at back of nebulizers

First Glass Concentric design that matches

Burgener or ESI’s PFA constant ID

More secure gas line attachments.

2012 Precision Glassblowing



**Optional Tygon tubing or Ezylok gas line fitting.
Wide range of models
Design a mix of Meinhard and GE features**

John A. Burgener


 Burgener Research Inc.

2004 Epond's Tyfoon Concentrics





the EPOND 'Tyfoon'

- *higher efficiency*
- *higher sensitivity*
- *higher salt tolerance*



**They have a Laval-type nozzle at the tip of the nebulizer,
Gives a more complete nebulization at higher sample flow rates.
Available in glass, quartz and PTFE.**

Epond's Lucida: 1 $\mu\text{l}/\text{min}$ Microconcentric



**The Lucida is unique
in that its natural aspirated uptake is below 10 $\mu\text{l}/\text{min}$.
It can run from 1 to 50 $\mu\text{l}/\text{min}$.**

John A. Burgener

 Burgener Research Inc.

Until 2001, all pneumatic nebulizers used induction.


Gas moving out of an orifice moves in a reasonably straight line. There is little sideways motion.

Molecules beside the gas stream are not pushed away, so they move towards the gas stream as though there was a vacuum.

As molecules move into the gas stream, they are impacted by the gas stream and move away with the gas stream.

Induction pulls liquid into gas stream.

Concentric
Cross flow
Babington V Groove
Original Parallel Path
All use induction



John A. Burgener

 Burgener Research Inc.

New concepts in nebulizers since 2000:

Burgener Enhanced Parallel Path Method

Evolved from Parallel Path

OneNeb Flow Blurring

Evolved from Concentric

Pizeoelectric Vibrating Mesh

Evolved from Ultra sonic

John A. Burgener

 Burgener Research Inc.

New concepts in nebulizers since 2000:

Burgener Enhanced Parallel Path Method

Evolved from Parallel Path

OneNeb Flow Blurring

Evolved from Concentric

Pizeoelectric Vibrating Mesh


Evolved from Ultra sonic

**ALL OF THESE ARE NEW CONCEPTS IN NEBULIZATION
NOT JUST ENHANCEMENTS OF PREVIOUS DESIGNS
AND NONE OF THESE USE INDUCTION**

John A. Burgener

 Burgener Research Inc.

Pizeoelectric Vibrating Mesh



A variation on Ultrasonics

The vibrating membrane has micro holes in it

The sample enters from the back and is pushed out the holes as the membrane vibrates.


Makes a fine mist with droplet size set by hole size.

Requires NO GAS FLOW

John A. Burgener

 Burgener Research Inc.

Pizeoelectric Vibrating Mesh




If the droplets are $< 5\mu\text{m}$ then they are too small to stick to the walls of the chamber. the chamber remains dry, and 90-100% of the sample goes to the torch. With no gas flow required, very low Argon flows are possible.

Pizeoelectric Vibrating Mesh

Sample In

To Torch



Sample Out

Advantages:

Droplet size set by hole size. All droplets same size.

If running $< 15 \mu\text{l}/\text{min}$ then

100% vaporizes and goes to the torch.

Very low flow Argon - same as Laser ablation.


John A. Burgener

 Burgener Research Inc.

Pizeoelectric Vibrating Mesh

Sample In

To Torch



Sample Out

Disadvantages:

If running $> 40 \mu\text{l}/\text{min}$, torch is overloaded, requiring use of desolvator.

Volume behind membrane $>$ sample aspirated.


Holes $< 5 \mu\text{m}$ are easy to plug.

Holes $< 5 \mu\text{m}$ are easy to salt.

John A. Burgener

 Burgener Research Inc.

Pizeoelectric Vibrating Mesh Summary



Likely for ICP/MS applications only.

Will be a niche usage for applications:

needing minimal sample sent to torch.

having lots of sample to send to nebulizer.

clean solutions only.

low salts only.

Main New nebulizers since 2000:

Burgener Enhanced Parallel Path Method

OneNeb Flow Blurring

Pizeoelectric Vibrating Mesh

**ALL OF THESE ARE NEW CONCEPTS IN NEBULIZATION
NOT JUST VARIATIONS ON PREVIOUS DESIGNS**

John A. Burgener

 *Burgener Research Inc.*

2001 Burgener Mira Mist




Looks like a Micro 3 but very different.

John A. Burgener


 Burgener Research Inc.

2001 Burgener Mira Mist

**Micro 3
Parallel Path Design**



**Mira Mist
Enhanced
Parallel Path Design**




**Dramatically different gas orifice
Dramatically different interaction
between gas and liquid**



Close-up of the Tip of a Burgener Nebulizer.

The Parallel Path Method

Induction pulls liquid into gas stream.




John A. Burgener

 Burgener Research Inc.

Enhanced Parallel Path Method

Surface tension
pulls liquid into spout.
Gas molecules impact liquid
breaking it into small droplets.

Interaction is at
center of gas stream,
Energy transfer is
10X greater
than from edge.



John A. Burgener

 Burgener Research Inc.

Enhanced Parallel Path Method

The Enhanced Method does not use gravity or induction. As such, it has zero back pressure and zero suction, and operates equally well in any orientation.


The liquid interacts with the gas stream in the central portion of the gas stream. Gas streams in capillaries have velocity gradients across the diameter of the capillary, with the slowest gas moving at the edges of the capillary and the fastest at the center. The central gas flow is 3 to 10 times as fast as the gas at the edge of the stream. If the central portion of the gas stream can impact the liquid, then the gas impacts the liquid with much more energy. Energy is related to the square of the velocity, so 3 to 10 times the speed is 9 to 100 times the energy. With such an increase in energy transfer from the gas to the liquid, the liquid is broken up into much smaller droplets. This produces a mist with average droplet sizes much smaller than any other method for the same gas flow and pressures.

John A. Burgener

 Burgener Research Inc.

Mira Mist Nebulizer Characteristics:

Normalized Size distribution for aerosol from Mira Mist on pure water, the most difficult liquid for a Mira Mist.




John A. Burgener

 Burgener Research Inc.

Data from Akbar Montaser,
George Washington University

Mira Mist Nebulizer Characteristics:

Normalized Size distribution for aerosol from Mira Mist on pure water, the most difficult liquid for a Mira Mist.



John A. Burgener


 Burgener Research Inc.

Data from Akbar Montaser,
George Washington University

High Stability

Mira Mist Stability Plot

< 1% RSD for 3% NaCl over a 4 hour period



John A. Burgener

 Burgener Research Inc.

Data from Jobin Yvon, France

**ENHANCED PARALLEL PATH
NEBULIZERS
for
ICP & ICP / MS
2012 Product line**



BETTER THAN THE BEST CONCENTRIC


**Excellent Sensitivity and Precision
with Little or No
Drifting, Salting, Plugging,
or Maintenance.**

BR *Burgener Research*

3 MONTH SATISFACTION GUARANTEE

Comparisons of sample passages


ARI MIST
225 μm




MIRA MIST
500 μm




T2100
760 μm




Enya Mist
60 μm



Mira Mist CE
500 μm



Glass Concentric
200 μm



Enhanced Parallel Path Summary



Ideal nebulizers for almost all uses:

Excellent stability and sensitivity.

Handle high to saturated salts.

Handle large particles.

Micro flow and macro flow.

Non contaminating fittings and materials.

Rugged: provides exceptionally long life.

John A. Burgener

 Burgener Research Inc.

Main New nebulizers since 2000:

Burgener Enhanced Parallel Path Method

OneNeb Flow Blurring

Pizeoelectric Vibrating Mesh


**ALL OF THESE ARE NEW CONCEPTS IN NEBULIZATION
NOT JUST VARIATIONS ON PREVIOUS DESIGNS**

John A. Burgener

 *Burgener Research Inc.*



OneNeb, a novel nebulizer for modern times



The new line of nebulizers

OneNeb is a new line of nebulizers based on the Flow Blurring® nebulization technology and specifically designed for liquid sample uptake, utilizing inductively-coupled plasma-based analytical techniques (e.g., optical emission spectrometry ICP-OES and mass spectrometry ICP-MS) for elemental chemical analysis.

Last News

Biomedal and Ingeniatics signs a distribution agreement for Cellena® bioencapsulator.

[Catálogo Interactivo de la OTRI](#)

Presentation

Our product

OneNeb is a new pneumatic nebulizer specially designed for use in liquid sample introduction systems for analytical atomic spectrometry techniques based on inductively coupled plasma (ICP-OES and ICP-MS).



The nebulization mechanism used is based on **Flow Blurring[®]** technology which is patented by the Spanish company Ingenialrics.

For Challenging Applications

– OneNeb Nebulizer

- Robust PFA and PEEK construction
 - Inert - resistant to strong acids such as HF
 - Resistant to breakage
 - Molded plastic design provides improved nebulizer to nebulizer reproducibility
- Constant diameter narrow bore tubing through to nebulizer tip
 - Ideal for high solids/particulates
 - Improved tolerance to high TDS samples
- Narrow aerosol size distribution provides improved precision
- Handles a wide flow range from 0.1 to 2 mL/min.
 - No sensitivity loss at low flow rates



Principle of Operation – Inert OneNeb Nebulizer

Tip geometry dimensioned to allow carrier gas to mix with the sample

Turbulent mixing of liquid and gas occurs



- Produces aerosol with a narrow size distribution range

Droplets mixed into gas flow

No area of low pressure



Unique Flow Blurring action increases nebulization efficiency

- Greater efficiency than conventional concentric nebulizer
- Improved sensitivity



Flow Blurring This is NOT using induction !


Concentric




John A. Burgener

 Burgener Research Inc.

Close Up Photos Flow Blurring



Glass Concentric




John A. Burgener

 Burgener Research Inc.

Concentric Photo courtesy of Geoff Coleman, MEINHARD®


OneNeb Normalized Stability Plot for 10 ppm S21+K solution over an 8 hour period



Agilent Application note:
“Analysis of wear metals and contaminants
in engine oils using the 4100 MP-AES”

Authors:
Phil Lowenstern and Elizabeth Reisman
Agilent Technologies, Melbourne, Australia

OneNeb Normalized Stability Plot for 10 ppm S21+K solution over an 8 hour period




Agilent Application note:
"Analysis of wear metals and contaminants
in engine oils using the 4100 MP-AES"

Authors:
Phil Lowenstern and Elizabeth Reisman
Agilent Technologies, Melbourne, Australia

OneNeb Normalized Stability Plot for 10 ppm S21+K solution over a 4 hour period


(This is half of the previous data, so 4 hours not 8.)



Agilent Application note:
“Analysis of wear metals and contaminants
in engine oils using the 4100 MP-AES”


Authors:
Phil Lowenstern and Elizabeth Reisman
Agilent Technologies, Melbourne, Australia

Mira Mist Normalized Stability Plot for 3% NaCl over a 4 hour period



DuraMist Long-term Stability

3.5%



OneNeb Review



New Concept, claims to have smaller droplets.

Pro:

High Precision, Low Detection Limits,

Because it has no low pressure zone, it should not salt.

Handles High Salts, Wide range of flow rates.

As with any micro nebulizer,

for < 50 ul/min can deliver nearly 100% to torch.

Con:

Small sample capillary - particles will plug

Rubber in Gas line fitting

- will cause some contaminations

Back Pressure can blow off sample lines

Smaller droplets will drown torches

unless running low flow or using a desolvator.


Claims of perfection seem to be based on

comparisons with non-optimized other nebulizers.



What's New in Nebulizers

PART 3: Some comments on Micro Nebulizers



The main Micro Nebulizers available today

Microconcentric Nebulizers:

- Aspire** – CETAC (TFE, PFA)
- CEI-100** – CETAC
- MicroMist** – Glass Expansion (glass)
- MicroFlow** – Elemental Scientific (PFA)
- HEN** – Meinhard (glass, quartz)
- Lucida, Zefyr, Boreal** – Epond (PTFE)

Non-Concentric Nebulizers:

- Mira Mist CE, Ari Mist, Enya Mist, PFA 250** – Burgener
- OneNeb** – Ingeniatics

Enya Mist: 0.2 - 50 $\mu\text{l}/\text{min}$ @ 45 psi
Ari Mist: 5 - 1500 $\mu\text{l}/\text{min}$ @ 45 psi
PFA 250: 50 - 1500 $\mu\text{l}/\text{min}$ @ 45 psi
Mira Mist CE: 1 - 2500 $\mu\text{l}/\text{min}$ @ 80 - 90 psi
Lucida: 1 - 40 $\mu\text{l}/\text{min}$ @ 35 psi
Zefyr: 5 - 100 $\mu\text{l}/\text{min}$ @ 35 psi

Boreal: 50 - 700 $\mu\text{l}/\text{min}$ @ 35 psi
OneNeb: 10 - 2000 $\mu\text{l}/\text{min}$ @ 40 psi
MicroMist: 100 - 500 $\mu\text{l}/\text{min}$ @ 30 - 60 psi
MicroFlow: 20 - 700 $\mu\text{l}/\text{min}$ @ 60 - 70 psi
High Efficiency Nebulizer: 30 - 300 $\mu\text{l}/\text{min}$ @ 90 - 170 psi
Microconcentric Nebulizer: 100 - 500 $\mu\text{l}/\text{min}$ @ 60 - 90 psi

Advantages of Micro and Nano flows

1 L of dry Argon can absorb 15 μl of water.

Running less than 15 $\mu\text{l}/\text{min}$ enables:

(conditional on droplet sizes being $< 10 \mu\text{m}$ and enough time allowed for vaporization or heating / desolvation being added)

ALL of the sample vaporizes and goes to torch.

The chamber stays dry

Washout times drop from minutes to milliseconds.

Torch stability improves.

Speciation effects dramatically reduced.

ICP/MS interferences drop dramatically.

MS cones/skimmers last much longer.


HPLC, μHPLC and CE can have ALL sample to torch.

Micro samples can be run over a long period of time.

Organic solvents can be run without effecting the plasma.

Transport Efficiency vs. Uptake

Chart assumes droplets have average sizes $< 10 \mu\text{m}$



John A. Burgener

 Burgener Research Inc.

Data from Geoff Coleman, Meinhard

Disadvantages of Micro Nebulizers?

Micro Concentrics usually have micro capillaries.

Micro particles can plug them.

Micro salts can plug them.

Micro damage can break them.

Advantages of Dry Samples

Torches can handle 40 - 60 $\mu\text{l}/\text{min}$ of sample

The variation probably relates to how much sample has vaporized.

**With Desolvation, Heating, Cooling,
some of the water can be removed.**

**With the sample all vaporized,
one gains the same benefits of low flow.**

**With High Efficiency Nebulizers, TOO MUCH sample
will arrive at the torch unless water is removed.**


More than 40-60 $\mu\text{l}/\text{min}$ water will quench / kill the plasma.

If you can dry 1 ml/min, the signal is 25X better than 40 $\mu\text{l}/\text{min}$.

Since 1988 Cetac Ultrasonic Nebulizers have used Desolvators

Efficiency approaches 30%, Sensitivity improves ~10x
Droplet size < 5 microns.
Potentially heavy solvent load so Desolvation essential.
Peltier cooling essential.
Membrane separator available and usually essential.
Desolvation interferences occur (eg., As III vs. As IV)

Does not handle high solids well.
Long washout times required due to desolvators.



Evaluation of an Ultrasonic Nebulizer Using Perkin-Elmer Sequential ICP Instrumentation - Robert Thomas, Cindy Anderau - Atomic Spectroscopy, 10 (2), 71 (1989). This paper describes the capabilities of the first commercially available ultrasonic nebulizer, the CETAC USN-5000 with two of Perkin Elmer's ICP optical emission spectrometers, the Plasma 40 and Plasma II. Various criteria were evaluated including detection limit, precision and memory test. The conclusion of the work was that the early design of this product showed great promise but required some improvements to be accepted and used as a truly routine analytical tool.

Desolvation Accessories

Now available for most nebulizers, not just Ultra Sonic

Sensitivity improves ~10x
Uses any 6 mm nebulizer
and cyclonic spray chamber

Desolvation interferences problematic
High solids still problematic

Available in HF-resistant versions
Available with membrane desolvation

Aridus II from CETAC
Apex from Elemental Scientific



John A. Burgener

 Burgener Research Inc.

Spray Chamber Temperature Control

Sensitivity improves
Uses any 6 mm nebulizer
and cyclonic spray chamber

Reduces solvent loading
Reduces oxide interferences in ICP/MS

Available in HF-resistant versions

PC3 from Elemental Scientific
IsoMist from Glass Expansion



John A. Burgener

 Burgener Research Inc.

Disadvantages of Micro Nebulizers?

Micro Concentrics usually have micro capillaries.

Micro particles can plug them.

Micro salts can plug them.


Micro damage can break them.

Micro Flow takes Macro time

to get from sample to nebulizer.

Flow Injection

Allows rapid sample turnaround for micro flows.



Low liquid flow rates mean very low sample throughput

Flow injection / fast rinse
Makes low flow practical.
Sample is quickly sucked up by vacuum, then slowly sent to nebulizer.

John A. Burgener

 Burgener Research Inc.

Available from CETAC, ESI, GE

Split Flow

Allows rapid sample turnaround for micro flows.

ENYA MIST - 0.2 to 50 $\mu\text{l}/\text{min}$
SPLIT FLOW CAPABILITY


Input 500 $\mu\text{l}/\text{min}$

Max flow of 50 $\mu\text{l}/\text{min}$
Desired flow 10 $\mu\text{l}/\text{min}$



Output to drain
490 $\mu\text{l}/\text{min}$


ENYA MIST - 0.2 to 50 $\mu\text{l}/\text{min}$




60 micron ID in nebulizer

SPLIT FLOW CAPABILITY

Similar configuration to
Mira Mist CE



**Enya Mist with
PE Elan 9000
ESI Cooled cyclonic**



Disadvantages of Micro Nebulizers?

Micro Concentrics usually have micro capillaries.

Micro particles can plug them.

Micro salts can plug them.

Micro damage can break them.

Micro Flow takes Macro time

to get from sample to nebulizer.

Long narrow capillaries cause signal broadening

when connected to LC / CE

LC, HPLC Connections

readily available from most manufacturers



Mira Mist LC interface



Meinhard and ESI standard fittings attach to LC




Glass Expansion LC interface



John A. Burgener

 Burgener Research Inc.

Signal Broadening Problem





John A. Burgener

 Burgener Research Inc.

CE Interfaces

Only available from Burgener and Cetac

Mira Mist CE interface




CE Capillary extends to tip: No signal broadening

John A. Burgener

 Burgener Research Inc.

Cetac CEI 100



Short length micro concentric: minimal signal broadening


Mira Mist CE: Interface for ICP/MS, CE & LC
Convenient Nano Sample Introduction
2 $\mu\text{l}/\text{min}$ to 2,500 $\mu\text{l}/\text{min}$




ICP/MS Interface for CE and LC
0.002 to 2.5 ml/min sample flow rates.
80+ psi for 1 L/min

Ideal Nebulizer for CE and LC interfacing to ICP/MS

The CE or LC capillary runs through the nebulizer body to the tip, enabling “instant” flow from the CE/LC capillary to the gas stream. There is no loss of signal as there is with nebulizers interfacing at the back of the nebulizer. The Mira Mist CE has zero suction or back pressure, so there is no interference with the flow from the CE or LC, just 100% transfer, instantly.




Detail of tip of Mira Mist CE




Mira Mist CE

*A workable Routine CE - ICP/MS interface
and Nano Sample Introduction Method*



Atomization is independent of sample flow.


As long as the make up solution covers the CE capillary, it maintains good electrical contact and the CE sample is washed directly to the gas orifice.



Close Up of Tip of Mira Mist CE



**Enhanced
Parallel Path
Method
with
CE
Capillary**




Other approach to μ HPLC interfacing

Signal broadening occurs when the sample passes from a tiny capillary into a wider capillary

Solution: Make a nebulizer with a similar size capillary.

The ENYA MIST: 60 μm ID

ENYA = Irish name meaning “Little Fire”.
The Enya Mist is designed to send
very little mist into the plasma (fire).



200 nl to 50 $\mu\text{l}/\text{min}$ Micro Flow Nebulizer
For nano HPLC interfacing
And for Split Flow applications

Standard sizes of micro-HPLC columns

(From G & T Septeck)

Column ID (mm)	Column Category	Flow Rate Range ($\mu\text{L}/\text{min}$)	Optimum Flow Rate ($\mu\text{L}/\text{min}$)	Analyte Capacity	Bed volume 15 cm column (μL)
1.0	microbore	20 - 200	35	$\sim 10 \mu\text{g}$	120
0.5	capillary	5 - 50	8	$\sim 2.5 \mu\text{g}$	30
0.3	capillary	2 - 20	3	$\sim 1 \mu\text{g}$	10
0.1	nano	0.25 - 2.5	0.5	$\sim 100 \text{ ng}$	2
0.075	nano	0.1 - 1	0.3	$\sim 25 \text{ ng}$	0.7

All micro and nano HPLC flow ranges covered by Enya Mist




Burgener Research Inc.




What's New in Nebulizers

PART 4: New 2013 nebulizers from Burgener The End is very near.



The newest Burgener Nebulizers (2013)

The X 175




Made in the same style as the Mira Mist CE
Greater stability at low flows than the Ari Mist

Wetted surfaces are Peek or PTFE
The body is strong PEEK.

Runs from 50 μ l/min to 1 ml/min
175 micron sample line ID.

The newest Burgener Nebulizers (2013)



The X series can be made in any size internal capillary,
from 60 - 760 micron ID

Which matches the main product line Enya Mist to T2100
But runs at higher gas pressures, lower gas flows, and greater stability.

They are expected to mainly be made
in the Ari Mist range of 100 - 250 micron ID.

The newest Burgener Nebulizers (2013)

The PFA 250



**All wetted surfaces are PFA.
The body is strong PEEK.**

**Runs from 50 μ l/min to 1.5 ml/min
250 micron sample line ID**

The newest Burgener Nebulizers (2013)

The PFA 250


1/16" PFA capillary
extends to the
tip of the nebulizer

All PFA in chamber



All wetted surfaces are PFA.
The body is strong PEEK.

Runs from 50 μ l/min to 1.5 ml/min
250 micron sample line ID



Thank you for your attention.